

# Medifil® II

## Collagen Particles



### Description

Medifil® II is 100% non-hydrolyzed type 1 bovine native collagen. Human BioSciences, Inc.'s proprietary Kollagen™ technology process protects and retains significantly more native triple helical protein structure, thus allowing superior stability of the molecule and scaffolding through all four phases of wound healing.

### Indications

For use and management of *burns, scrapes, blisters, sores, ulcers, acute and chronic wounds, superficial, partial and full thickness wounds, infected and non-infected wounds, and minimal to heavily exudating wounds*. Easy application and product compatibility with topical agents and or dressings.

### Application

Cleanse the wound per physician order. Apply Medifil® II particles ¼ inch to entire wound bed surface, including any undermining and or tunneling. Cover the wound with a secondary absorbent dressing to promote thermoregulation and maintain an optimal moist wound healing environment. Frequency of dressing change is based on exudate amount, clinical findings, and per physician recommendations.

### Contraindications

- ▲ Known bovine sensitivities.
- ▲ Third-degree burns



- Single-patient use
- Sterile until opened
- Store at room temperature
- Do not use if package is damaged

Manufactured in the USA  
Human BioSciences, Inc.  
940 Clopper Rd. Gaithersburg, MD 20878 USA  
[info@humanbiosciences.com](mailto:info@humanbiosciences.com)  
+1 301-740-1893  
©2020

